AGS

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Contents

For Your Safety and that of Your Patients.....	18
Intended Use.....	19
System description	19
Preparation.....	20
for anaesthetic workstations	20
for intensive care ventilators	23
Operation.....	24
Shut-down.....	24
Care.....	25
Stripping down.....	25
Cleaning/disinfecting	25
Sterilizing.....	25
Maintenance	26
Replacing particle filter	26
Replacing flow indicator.....	26
Fault, Cause, Remedy	27
What's What	28
Technical Data	29
Order List.....	30
Order List of the Anaesthetic Gas Scavenging System AGS.....	31

For Your Safety and that of Your Patients

Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observation of these instructions. The apparatus is only to be used for purposes specified here.

Maintenance

The apparatus must be inspected and serviced regularly by trained service personnel at six monthly intervals (and a record kept).

Repair and general overhaul of the apparatus may only be carried out by trained service personnel.

We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them. Only authentic Dräger spare parts may be used for maintenance.

Observe chapter "Maintenance Intervals".

Accessories

Do not use accessory parts other than those in the order list.

Safety inspections

The apparatus must be submitted to regular safety inspections.

Liability for proper function or damage

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus is serviced or repaired by personnel not employed or authorized by DrägerService or if the apparatus is used in a manner not conforming to its intended use.

Dräger cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Dräger are likewise not modified by the recommendations given above.

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Intended use

AGS - Anaesthetic Gas Scavenging System

For extracting and removing excess anaesthetic gas and expiratory gas.

For reducing the concentration of anaesthetic agent around anaesthetic workstations.

For scavenging contaminated expiratory gases during ventilation.

System description

The anaesthetic gas scavenging system consists of:

- Transfer system: corrugated hoses which connect the exhaust ports on the anaesthetic machine/ventilator to the receiving system.
- Receiving system: head with buffer container, built-in flow indicator and scavenging hose.
- Disposal system: terminal unit, ejector and exhaust air line.

The Dräger anaesthetic gas scavenging system operates as an "open system". The excess anaesthetic gases/expiration gases are channelled into the buffer container of the receiving system, from where they are exhausted in a continuous stream into the disposal system.

The openings in the receiving system prevent any increase in expiratory pressure or any negative pressure in the breathing system.


Preparation

for anaesthetic workstations

Assembly: mounting AGS


Using the built-in bracket on the anaesthetic system
example Cicero EM

- Hook the receiving system on to the bracket on the back panel:
Fit the slots of the receiving system over the two holders and allow the receiving system to slide down into place.


Assembly on anaesthetic machine trolley
example Sulla

- Use appropriate AGS bracket for anaesthetic system.
- Hang receiving system on the bracket on trolley.

Assembly on rail


- Fit the rail clamp to the receiving system.
- Hook the receiving system to a rail and secure with the twist-knob.

Connecting AGS

Connecting the scavenging hose


- 1 Fit the scavenging (exhaust) hose to the port on the receiving system.
- 2 Connect the anaesthetic exhaust connector to the scavenging hose.

Connecting transfer hoses


Anaesthetic machines with **one** exhaust port (waste-gas socket), example Cicero EM

- 3 Route transfer hose from the rear through the opening on anaesthetic machine and connect hose to waste-gas socket on breathing system.
- 4 Connect transfer hose to socket on receiving system.
- Make sure that the connections are not leaking.
- 5 Make sure that second transfer hose socket is closed with a screw cap.

Anaesthetic machines with **two** exhaust ports
(waste-gas sockets), example Sulla 909 with Ventilog


- Unscrew cap from second transfer hose socket on the receiving system.
- Connect both transfer hoses.

Connecting the sample-gas return hose

When the extracted sample gas is not returned,
example PM 8050 and Circle System 8


- 1 Connect hose to outlet for measured gas on monitor.
- 2 Plug the hose connector into the socket on the receiving system until it clicks into place.

for intensive care ventilators
Evita 4, Evita 2 dura

AGS Assembly


- Fit the rail clamp to the receiving system.
- 1 Hook the receiving system to a rail and secure with the twist-knob.

Connecting AGS


Connecting the transfer hose

- 2 Connect the transfer hose to the flow sensor of the ventilator by means of the adapter.
 - 3 Insert the transfer hose in the port of the receiving system.
- Make sure that the port for the second transfer hose is locked by a sealing screw.


Operation

- 1 Push probe into terminal unit until it engages.
- 2 The indicator will turn green and scavenging will commence.


- 3 The float in the flow indicator should stay between the upper and lower marks. If it does not, see section on "Fault, Cause, Remedy", page 27.

Do not block the openings in the receiving system; otherwise, a shortage of fresh gas may be caused in the breathing system.

Shut-down

- Pull black casing on probe back to disengage it. Scavenging will be switched off.
- The indicator will no longer be green.

Care

Stripping down

For cleaning purposes only

- Disconnect transfer hoses from receiving system and from anaesthetic machine.
- Disconnect scavenging hose.
- Disconnect the sample-gas return hose, if used. Pull back the ring of the coupling, and the plug will be ejected.
- Disconnect the receiving system from the anaesthetic machine/ventilator.
- Rotate buffer container by 45° and remove.
- Remove the flow indicator and filter if the receiving system head is to be cleaned and disinfected in the autoclave.

Recommended cleaning routines:

- transfer hoses: weekly
- scavenging hose: monthly
- buffer container: weekly
- receiving system head: monthly

Cleaning/disinfecting

Use surface disinfectants for disinfection.

For reasons of material compatibility, use disinfectants based on

- aldehydes
- alcohols
- quaternary ammonium compounds.

Do not use

- phenols
- halogen-releasing compounds
- strong organic acids
- oxygen-releasing compounds
- ethylene oxide sterilization.

We recommend that disinfectants on the current DGHM (DGHM: German Society for Hygiene and Microbiology) list are used. The DGHM list (mhp-Verlag GmbH, Wiesbaden) gives the composition of each disinfectant. For countries where the DGHM list is not familiar, we recommend the types of disinfectants given above.

Wipe disinfection

- Wipe off any coarse impurities with a disposable cloth.

Receiving system head, flow indicator, terminal unit, sample-gas return socket:

- Wipe disinfect, with Buraton 10 F or Terralin (Schülke & Mayr, Norderstedt), for example. Observe manufacturer's instructions.
- Do not allow any liquid to get into terminal units.

Clean and disinfect in cleaning and disinfecting machine

Transfer hoses, scavenging hose, buffer container, receiving system head (without flow indicator and filter):

- Disinfect with moist heat up to no more than 95 °C, for at least 10 minutes. Do not add any alkaline or chlorine-releasing cleaning and disinfecting agents.

Sterilizing

Carry out when required.

Buffer container

- sterilize in steam up to 134 °C.

Before re-using

- Re-assemble, see page 20 to page 24.

Test of functioning

- Plug the probe into the terminal unit: the float in the flow indicator should lie between the upper and lower marks. If it does not, see section on "Fault, Cause, Remedy", page 27.


Maintenance

Clean and disinfect equipment and parts before any maintenance procedures – and before returning for repair.

Particle filter	replace if blocked
Flow indicator	replace if float sticks
Inspection and service	every six months by trained service personnel
Safety inspection	every 6 months


Replacing particle filter

- Remove buffer container.
- 1 Loosen retaining nut and remove blocked particle filter from the flow indicator.
 - 2 Replace with new particle filter and re attach retaining nut.
- Dispose of used particle filter with normal waste.
 - Put buffer container back.

Replacing flow indicator

- Remove buffer container.
- 3 Loosen retaining nut.
 - 4 Unscrew faulty flow indicator and
- fit new flow indicator.
 - Put container for buffer volume back.


Fault, Cause, Remedy

Fault	Cause	Remedy
Float sinks below bottom mark	Particle filter blocked.	Replace particle filter, page 26.
	Insufficient suction power of the vacuum generator.	Ask Dräger Service to set the correct vacuum generator suction for the AGS working range, or ask Dräger Service to replace the vacuum generator with a more powerful generator (conforming to EN737).
Float rises above top mark	No particle filter.	Install particle filter, page 26.
	Suction power of the vacuum generator too high.	Ask Dräger Service to set the correct vacuum generator suction for the AGS working range.
Float is stuck in flow indicator, position does not change during start-up	Faulty flow indicator.	Replace flowmeter, page 26.

What's What

Receiving system

- 1 Receiving system head
- 2 Socket for scavenging hose
- 3 Flow indicator with float to show scavenging flow
- 4 Particle filter
- 5 Buffer container
- 6 Socket for transfer hose
- 7 Socket for second transfer hose, closed with screw cap
- 8 Socket for sample-gas return hose

Technical Data

Ambient conditions

During operation:

Temperature	15 to 35 °C
Atmospheric pressure	700 to 1060 hPa
Rel. humidity	0 to 98 %, without condensation

During storage:

Temperature	–20 to 70 °C
Atmospheric pressure	500 to 1100 hPa
Rel. humidity	0 to 98 %, without condensation

Performance data

Effect on the anaesthetic system/ventilator:

Increase in exhalation resistance	<0.5 mbar at 30 L/min <3.5 mbar at 75 L/min
Generated vacuum pressure	<1 Pa (0.01 mbar)
Generated flow	<50 mL/min

Effect on the environment:

Gas loss, dependent on fresh gas flow and ventilation parameters	<25 mL/min (typical)
--	----------------------

Operating data

Operating capacity of scavenging flow
Setting with service flow indicator

30 to 50 L/min
65 to 75 L/min

Particle filter pore width
Particle filter resistance

0.1 mm
5 mbar at 50 L/min

Supply performance requirements
across a flow resistance of 2 kPa
across a flow resistance of 1 kPa

min. 25 L/min suction flow
max. 50 L/min suction flow

Dimensions (W x H x D) mm

120 x 350 x 200

Weight, without hoses

ca. 970 g

Order List

Name	Order No.
AGS receiving system	M 33 300
AGS transfer hose 0,5 m	M 35 016
AGS transfer hose 1 m	M 33 295
AGS transfer hose 1,5 m	M 33 296
AGS scavenging hose 0,5 m	M 35 015
AGS scavenging hose 1,5 m	M 33 297
AGS scavenging hose 3 m	M 33 298
AGS scavenging hose 5 m	M 33 299
Anaesthetic waste gas probe	G 60 440
Clamp modification kit	M 32 967
AGS bracket for Cicero / Sulla	M 33 303
AGS bracket for Titus / SA 2	M 32 976
AGS adapter for ventilators	84 13 433
Hose set for sample-gas return	M 32 692
Spare parts	
Buffer container	M 33 292
Flow indicator	M 33 293
Screw cap	M 33 291
Particle filter	M 33 294

Order List of the Anaesthetic Gas Scavenging System AGS

	AGS- receiving system	AGS transfer hose 0.5 m	AGS- transfer hose 0.75 m	AGS- transfer hose 1.0 m	AGS- transfer hose 1.5 m	AGS- scavenging hose 0.5 m	AGS- scavenging hose 1.5 m	AGS- scavenging hose 3.0 m	AGS- scavenging hose 5.0 m	Anaes- thetic waste gas probe	AGS- adapter	Clamp modifica- tion kit	AGS- bracket for Cicero / Sulla M 33 303	AGS- bracket for Titus/ SA 2 M 32 976
	M 33 300	M 35 016	M 33 955	M 33 295	M 33 296	M 35 015	M 33 297	M 33 298	M 33 299	G 60 440	84 13 433	M 32 967	M 33 303	M 32 976
Cato	X				X		A	A	A	X				
Cato ceiling-mounted	X				X		X			X				
Cicero	X			T	T		A	A	A	X			X	
Cicero ceiling-mounted	X			T	T		X			X			X	
Cicero EM	X			X			A	A	A	X				
Cicero EM ceiling-mounted	X			X			X			X				
Julian	X		X				A	A	A	X				
Julian wall-mount. (new installation)	X	X				X				X				
Julian wall-mount. (Replacement)	X	X					A	A	A	X				
Julian ceiling-mounted	X	X					A	A	A					
Livius	X			X (2x)			A	A	A	X		X		
Remus	X			X	X		A	A	A	X		X		
SA 2 without absorber ¹⁾	X			X (2x)			A	A	A	X		R		R
SA 2 with absorber ¹⁾	X			X			A	A	A	X		R		R
Sulla	X			X (2x)			A	A	A	X		R	R	
Titus	X			X			A	A	A	X		X		
Titus with trolley	X			X (2x)			A	A	A	X		R		R
Trajan	X			X	X		A	A	A	X		X		
Evita 4	X		X				A	A	A	X	X	X		
Evita 2 dura	X		X				A	A	A	X	X	X		

Legend: X = absolute necessary
T = either transfer hose
A = either scavenging hose
R = either modification kit
1 = only for SA 2 additional necessary: adapter M 29 430

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Directive 93/42/EEC
concerning Medical Devices

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☞ Moislinger Allee 53 – 55

D-23542 Lübeck

☎ +49 451 8 82 - 0

☎ 26 80 70

FAX +49 451 8 82-20 80

🌐 <http://www.draeger.com>

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Germany

☞ Moislinger Allee 53 – 55

D-23542 Lübeck

☎ +49 451 8 82 - 0

☎ 26 80 70

FAX +49 451 8 82-20 80

🌐 <http://www.draeger.com>

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Dräger Medizintechnik GmbH